


CORIN



DC-Fit™ Acetabular Cup System

UHMWPE inserts dialable to any one of twelve positions

Inserts available in standard or 10° offset configurations
with 32mm, 28mm, 26mm and 22mm inside diameters

Bi-Coated implants (plasma sprayed titanium and
hydroxyapatite) for cementless fixation

DC-FIT™ Ordering Information

Cementless Implants Bi Coated (with holes, no apical pin)

S193.142B	DC-Fit Bi Coated 42mm 2 holes (Size 1 Liner)
S193.144B	DC-Fit Bi Coated 44mm 2 holes (Size 1 Liner)
S193.246B	DC-Fit Bi Coated 46mm 2 holes (Size 2 Liner)
S193.248B	DC-Fit Bi Coated 48mm 2 holes (Size 2 Liner)
S193.350B	DC-Fit Bi Coated 50mm 3 holes (Size 3 Liner)
S193.352B	DC-Fit Bi Coated 52mm 3 holes (Size 3 Liner)
S193.454B	DC-Fit Bi Coated 54mm 3 holes (Size 4 Liner)
S193.456B	DC-Fit Bi Coated 56mm 3 holes (Size 4 Liner)
S193.558B	DC-Fit Bi Coated 58mm 3 holes (Size 5 Liner)
S193.560B	DC-Fit Bi Coated 60mm 3 holes (Size 5 Liner)
S193.662B	DC-Fit Bi Coated 62mm 3 holes (Size 6 Liner)
S193.664B	DC-Fit Bi Coated 64mm 3 holes (Size 6 Liner)
S193.666B	DC-Fit Bi Coated 66mm 3 holes (Size 6 Liner)
S193.668B	DC-Fit Bi Coated 68mm 3 holes (Size 6 Liner)
S193.670B	DC-Fit Bi Coated 70mm 3 holes (Size 6 Liner)

Cementless Implants Bi Coated (no holes, no apical pin)

S194.142B	DC-Fit Bi Coated 42mm (Size 1 Liner)
S194.144B	DC-Fit Bi Coated 44mm (Size 1 Liner)
S194.246B	DC-Fit Bi Coated 46mm (Size 2 Liner)
S194.248B	DC-Fit Bi Coated 48mm (Size 2 Liner)
S194.350B	DC-Fit Bi Coated 50mm (Size 3 Liner)
S194.352B	DC-Fit Bi Coated 52mm (Size 3 Liner)
S194.454B	DC-Fit Bi Coated 54mm (Size 4 Liner)
S194.456B	DC-Fit Bi Coated 56mm (Size 4 Liner)
S194.558B	DC-Fit Bi Coated 58mm (Size 5 Liner)
S194.560B	DC-Fit Bi Coated 60mm (Size 5 Liner)
S194.662B	DC-Fit Bi Coated 62mm (Size 6 Liner)
S194.664B	DC-Fit Bi Coated 64mm (Size 6 Liner)
S194.666B	DC-Fit Bi Coated 66mm (Size 6 Liner)
S194.668B	DC-Fit Bi Coated 68mm (Size 6 Liner)
S194.670B	DC-Fit Bi Coated 70mm (Size 6 Liner)

Cementless Implants Bi Coated (with holes and apical pin)

148.142B	DC-Fit Bi Coated 42mm 4 holes (Size 1 Liner)
148.144B	DC-Fit Bi Coated 44mm 4 holes (Size 1 Liner)
148.246B	DC-Fit Bi Coated 46mm 4 holes (Size 2 Liner)
148.248B	DC-Fit Bi Coated 48mm 4 holes (Size 2 Liner)
148.350B	DC-Fit Bi Coated 50mm 12 holes (Size 3 Liner)
148.352B	DC-Fit Bi Coated 52mm 12 holes (Size 3 Liner)
148.454B	DC-Fit Bi Coated 54mm 12 holes (Size 4 Liner)
148.456B	DC-Fit Bi Coated 56mm 12 holes (Size 4 Liner)
148.558B	DC-Fit Bi Coated 58mm 12 holes (Size 5 Liner)
148.560B	DC-Fit Bi Coated 60mm 12 holes (Size 5 Liner)
148.662B	DC-Fit Bi Coated 62mm 12 holes (Size 6 Liner)
148.664B	DC-Fit Bi Coated 64mm 12 holes (Size 6 Liner)
148.666B	DC-Fit Bi Coated 66mm 12 holes (Size 6 Liner)
148.668B	DC-Fit Bi Coated 68mm 12 holes (Size 6 Liner)
148.670B	DC-Fit Bi Coated 70mm 12 holes (Size 6 Liner)

Cementless Implants Bi Coated (with holes and apical pin)

190.142B	DC-Fit Bi Coated 42mm 2 holes (Size 1 Liner)
190.144B	DC-Fit Bi Coated 44mm 2 holes (Size 1 Liner)
190.246B	DC-Fit Bi Coated 46mm 2 holes (Size 2 Liner)
190.248B	DC-Fit Bi Coated 48mm 2 holes (Size 2 Liner)
190.350B	DC-Fit Bi Coated 50mm 3 holes (Size 3 Liner)
190.352B	DC-Fit Bi Coated 52mm 3 holes (Size 3 Liner)
190.454B	DC-Fit Bi Coated 54mm 3 holes (Size 4 Liner)
190.456B	DC-Fit Bi Coated 56mm 3 holes (Size 4 Liner)
190.558B	DC-Fit Bi Coated 58mm 3 holes (Size 5 Liner)
190.560B	DC-Fit Bi Coated 60mm 3 holes (Size 5 Liner)
190.662B	DC-Fit Bi Coated 62mm 3 holes (Size 6 Liner)
190.664B	DC-Fit Bi Coated 64mm 3 holes (Size 6 Liner)
190.666B	DC-Fit Bi Coated 66mm 3 holes (Size 6 Liner)
190.668B	DC-Fit Bi Coated 68mm 3 holes (Size 6 Liner)
190.670B	DC-Fit Bi Coated 70mm 3 holes (Size 6 Liner)

UHMWPE Liners - Standard

150.122	Dialable Standard Liner 22mm Size 1
150.222	Dialable Standard Liner 22mm Size 2
150.226	Dialable Standard Liner 26mm Size 2
150.228	Dialable Standard Liner 28mm Size 2
150.322	Dialable Standard Liner 22mm Size 3
150.326	Dialable Standard Liner 26mm Size 3
150.328	Dialable Standard Liner 28mm Size 3
150.422	Dialable Standard Liner 22mm Size 4
150.426	Dialable Standard Liner 26mm Size 4
150.428	Dialable Standard Liner 28mm Size 4
150.432	Dialable Standard Liner 32mm Size 4
150.522	Dialable Standard Liner 22mm Size 5
150.526	Dialable Standard Liner 26mm Size 5
150.528	Dialable Standard Liner 28mm Size 5
150.532	Dialable Standard Liner 32mm Size 5
150.622	Dialable Standard Liner 22mm Size 6
150.626	Dialable Standard Liner 26mm Size 6
150.628	Dialable Standard Liner 28mm Size 6
150.632	Dialable Standard Liner 32mm Size 6

Cementless Implants Ti Porous Coated (with holes, no apical pin)

S193.142T	DC-Fit Ti Porous Coated 42mm 2 holes (Size 1 Liner)
S193.144T	DC-Fit Ti Porous Coated 44mm 2 holes (Size 1 Liner)
S193.246T	DC-Fit Ti Porous Coated 46mm 2 holes (Size 2 Liner)
S193.248T	DC-Fit Ti Porous Coated 48mm 2 holes (Size 2 Liner)
S193.350T	DC-Fit Ti Porous Coated 50mm 3 holes (Size 3 Liner)
S193.352T	DC-Fit Ti Porous Coated 52mm 3 holes (Size 3 Liner)
S193.454T	DC-Fit Ti Porous Coated 54mm 3 holes (Size 4 Liner)
S193.456T	DC-Fit Ti Porous Coated 56mm 3 holes (Size 4 Liner)
S193.558T	DC-Fit Ti Porous Coated 58mm 3 holes (Size 5 Liner)
S193.560T	DC-Fit Ti Porous Coated 60mm 3 holes (Size 5 Liner)
S193.662T	DC-Fit Ti Porous Coated 62mm 3 holes (Size 6 Liner)
S193.664T	DC-Fit Ti Porous Coated 64mm 3 holes (Size 6 Liner)
S193.666T	DC-Fit Ti Porous Coated 66mm 3 holes (Size 6 Liner)
S193.668T	DC-Fit Ti Porous Coated 68mm 3 holes (Size 6 Liner)
S193.670T	DC-Fit Ti Porous Coated 70mm 3 holes (Size 6 Liner)

Cementless Implants Ti Porous Coated (no holes, no apical pin)

S194.142T	DC-Fit Ti Porous Coated 42mm (Size 1 Liner)
S194.144T	DC-Fit Ti Porous Coated 44mm (Size 1 Liner)
S194.246T	DC-Fit Ti Porous Coated 46mm (Size 2 Liner)
S194.248T	DC-Fit Ti Porous Coated 48mm (Size 2 Liner)
S194.350T	DC-Fit Ti Porous Coated 50mm (Size 3 Liner)
S194.352T	DC-Fit Ti Porous Coated 52mm (Size 3 Liner)
S194.454T	DC-Fit Ti Porous Coated 54mm (Size 4 Liner)
S194.456T	DC-Fit Ti Porous Coated 56mm (Size 4 Liner)
S194.558T	DC-Fit Ti Porous Coated 58mm (Size 5 Liner)
S194.560T	DC-Fit Ti Porous Coated 60mm (Size 5 Liner)
S194.662T	DC-Fit Ti Porous Coated 62mm (Size 6 Liner)
S194.664T	DC-Fit Ti Porous Coated 64mm (Size 6 Liner)
S194.666T	DC-Fit Ti Porous Coated 66mm (Size 6 Liner)
S194.668T	DC-Fit Ti Porous Coated 68mm (Size 6 Liner)
S194.670T	DC-Fit Ti Porous Coated 70mm (Size 6 Liner)

Cementless Implants Ti Porous Coated (with holes and apical pin)

148.142T	DC-Fit Ti Porous Coated 42mm 4 holes (Size 1 Liner)
148.144T	DC-Fit Ti Porous Coated 44mm 4 holes (Size 1 Liner)
148.246T	DC-Fit Ti Porous Coated 46mm 4 holes (Size 2 Liner)
148.248T	DC-Fit Ti Porous Coated 48mm 4 holes (Size 2 Liner)
148.350T	DC-Fit Ti Porous Coated 50mm 12 holes (Size 3 Liner)
148.352T	DC-Fit Ti Porous Coated 52mm 12 holes (Size 3 Liner)
148.454T	DC-Fit Ti Porous Coated 54mm 12 holes (Size 4 Liner)
148.456T	DC-Fit Ti Porous Coated 56mm 12 holes (Size 4 Liner)
148.558T	DC-Fit Ti Porous Coated 58mm 12 holes (Size 5 Liner)
148.560T	DC-Fit Ti Porous Coated 60mm 12 holes (Size 5 Liner)
148.662T	DC-Fit Ti Porous Coated 62mm 12 holes (Size 6 Liner)
148.664T	DC-Fit Ti Porous Coated 64mm 12 holes (Size 6 Liner)
148.666T	DC-Fit Ti Porous Coated 66mm 12 holes (Size 6 Liner)
148.668T	DC-Fit Ti Porous Coated 68mm 12 holes (Size 6 Liner)
148.670T	DC-Fit Ti Porous Coated 70mm 12 holes (Size 6 Liner)

Cementless Implants Ti Porous Coated (with holes and apical pin)

190.142T	DC-Fit Ti Porous Coated 42mm 2 holes (Size 1 Liner)
190.144T	DC-Fit Ti Porous Coated 44mm 2 holes (Size 1 Liner)
190.246T	DC-Fit Ti Porous Coated 46mm 2 holes (Size 2 Liner)
190.248T	DC-Fit Ti Porous Coated 48mm 2 holes (Size 2 Liner)
190.350T	DC-Fit Ti Porous Coated 50mm 3 holes (Size 3 Liner)
190.352T	DC-Fit Ti Porous Coated 52mm 3 holes (Size 3 Liner)
190.454T	DC-Fit Ti Porous Coated 54mm 3 holes (Size 4 Liner)
190.456T	DC-Fit Ti Porous Coated 56mm 3 holes (Size 4 Liner)
190.558T	DC-Fit Ti Porous Coated 58mm 3 holes (Size 5 Liner)
190.560T	DC-Fit Ti Porous Coated 60mm 3 holes (Size 5 Liner)
190.662T	DC-Fit Ti Porous Coated 62mm 3 holes (Size 6 Liner)
190.664T	DC-Fit Ti Porous Coated 64mm 3 holes (Size 6 Liner)
190.666T	DC-Fit Ti Porous Coated 66mm 3 holes (Size 6 Liner)
190.668T	DC-Fit Ti Porous Coated 68mm 3 holes (Size 6 Liner)
190.670T	DC-Fit Ti Porous Coated 70mm 3 holes (Size 6 Liner)

UHMWPE Liners - 10° Offset

152.122	Dialable 10° Offset Liner 22mm Size 1
152.222	Dialable 10° Offset Liner 22mm Size 2
152.226	Dialable 10° Offset Liner 26mm Size 2
152.228	Dialable 10° Offset Liner 28mm Size 2
152.322	Dialable 10° Offset Liner 22mm Size 3
152.326	Dialable 10° Offset Liner 26mm Size 3
152.328	Dialable 10° Offset Liner 28mm Size 3
152.422	Dialable 10° Offset Liner 22mm Size 4
152.426	Dialable 10° Offset Liner 26mm Size 4
152.428	Dialable 10° Offset Liner 28mm Size 4
152.432	Dialable 10° Offset Liner 32mm Size 4
152.522	Dialable 10° Offset Liner 22mm Size 5
152.526	Dialable 10° Offset Liner 26mm Size 5
152.528	Dialable 10° Offset Liner 28mm Size 5
152.532	Dialable 10° Offset Liner 32mm Size 5
152.622	Dialable 10° Offset Liner 22mm Size 6
152.626	Dialable 10° Offset Liner 26mm Size 6
152.628	Dialable 10° Offset Liner 28mm Size 6
152.632	Dialable 10° Offset Liner 32mm Size 6



Corin
The Corinium Centre
Cirencester, GL7 1YJ
Tel: +44 (0) 1285 659 866
Fax: +44 (0) 1285 658 960
Email: corin@corin.co.uk
www.coringroup.com



Corin Japan KK
2-34-1-501 Higashi-Mikuni
Yodogawa-Ku
Osaka 532-0002, Japan
Tel: +81 6 6391 8651
Fax: +81 6 6391 8653
Email: welcome@corin-j.com



Corin USA
10500 University Center Drive
Suite 190, Tampa
Florida 33612, USA
Tel: +1 813 977 4469
Fax: +1 813 979 0042
Email: sales@corinusa.com



Alphanorm GmbH
Fischbacher Straße 92
D-66287 Quierschied, Germany
Tel: +49 (0) 6897 600 800
Fax: +49 (0) 6897 600 8050
Email: info@alphanorm.de
www.alphanorm.de